



HDX mouse mAb

Catalog No	YP-mAb-10138
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HDX CXorf43
Immunogen	Synthesized peptide derived from human HDX AA range: 321-371
Specificity	This antibody detects endogenous levels of human HDX
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Deltex1;E3 ubiquitin-protein ligase DTX1;EC:2.3.2.27;hDTX1;hDx 1;HDX;HDX CXorf43;Highly divergent homeobox
Observed Band	67 kDa
Calculated Molecular Weight	620 aa, 67 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images