



ZNH12 mouse mAb

Catalog No	YP-mAb-10142
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNHIT2 C11orf5
Immunogen	Synthesized peptide derived from human ZNH12 AA range: 26-76
Specificity	This antibody detects endogenous levels of human ZNH12
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZNH12;ZNHIT2 C11orf5;Zinc finger HIT domain-containing protein 2 (Protein FON)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Low expression in most tissues; highly expressed in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images