



UDB11 mouse mAb

Catalog No	YP-mAb-10149
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UGT2B11
Immunogen	Synthesized peptide derived from human UDB11 AA range: 53-103
Specificity	This antibody detects endogenous levels of human UDB11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UDB11;UGT2B11;UDP-glucuronosyltransferase 2B11 (UDPGT 2B11) (EC 2.4.1.17)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Microsome membrane ; Single-pass membrane protein . Endoplasmic reticulum membrane ; Single-pass membrane protein .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images