



LRGUK mouse mAb

Catalog No	YP-mAb-10152
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LRGUK
Immunogen	Synthesized peptide derived from human LRGUK AA range: 688-738
Specificity	This antibody detects endogenous levels of human LRGUK
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ32786;Leucine-rich repeat and guanylate kinase domain-containing protein; LRGUK;LRGUK_HUMAN
Observed Band	94kD
Calculated Molecular Weight	94kD
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome . Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, cilium basal body . Localizes to the acrosome and acroplaxome in round spermatids. Localizes to the manchette during spermiogenesis. Also found in the basal body of elongating spermatids, and in primary cilia of somatic cells. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images