



ZHX1 mouse mAb

Catalog No	YP-mAb-10155
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZHX1
Immunogen	Synthesized peptide derived from human ZHX1 AA range: 736-786
Specificity	This antibody detects endogenous levels of human ZHX1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	zinc fingers and homeoboxes 1;ZHX1;Zinc fingers and homeoboxes protein 1
Observed Band	100-120 kDa
Calculated Molecular Weight	873 aa, 98 kDa
Cell Pathway	Nucleus . Colocalized in the nucleus with DNMT3B.
Tissue Specificity	Ubiquitously expressed. Expressed in podocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images