



CS057 mouse mAb

Catalog No	YP-mAb-10156
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C19orf57
Immunogen	Synthesized peptide derived from human CS057 AA range: 246-296
Specificity	This antibody detects endogenous levels of human CS057
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CS057;C19orf57;Uncharacterized protein C19orf57 (Pre-T/NK cell-associated protein 3B3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Chromosome . During meiosis, recruited to chromosomes and localizes on recombination sites in a double-strand break-dependent manner. First appears on the chromosome axis at leptotene. Along with the progression of meiotic recombination, released from the axis to form bridge-like structures linking homolog axes before they are synapsed. Finally, located between synapsed homolog axes and on the synaptonemal complex (SC). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images