



MMSA mouse mAb

Catalog No	YP-mAb-10162
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ALDH6A1 MMSDH
Immunogen	Synthesized peptide derived from human MMSA AA range: 220-270
Specificity	This antibody detects endogenous levels of human MMSA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Antigen MMSA 1;Asp His His Cys domain containing protein 9;CGI89;CXorf11; DHHC-9;DHHC9;MMSA1;MRXSZ;Palmitoyltransferase ZDHHC9; ZDHC9_HUMAN;ZDHHC 9;ZDHHC10;ZDHHC9;Zinc finger DHHC domain containing 10;Zinc finger DHHC domain containing protein 9;Zinc finger DHHC domain-containing protein 9;Zinc finger DHHC type containing 9;Zinc finger protein 379;Zinc finger protein 380;ZNF379;ZNF380;MMSA;ALDH6A1 MMSDH; Methylmalonate-semialdehyde dehydrogenase [acylating];mitochondrial (MMSDH) (Malonate-semialdehyde dehydrogenase [acylating]) (EC 1.2.1.18) (EC 1.2.1.27) (Aldehyde dehydrogenase family 6 member A1); mitochondrial (MMSDH) (Malonate-semialdehyde dehydrogenase [acylating]) (EC 1.2.1.18) (EC 1.2.1.27) (Aldehyde dehydrogenase family 6 member A1)
Observed Band	41kD
Calculated Molecular Weight	41kDa
Cell Pathway	Mitochondrion.
Tissue Specificity	
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images