



MPCP mouse mAb

Catalog No	YP-mAb-10163
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC25A3 PHC OK/SW-cl.48
Immunogen	Synthesized peptide derived from human MPCP AA range: 102-152
Specificity	This antibody detects endogenous levels of human MPCP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MPCP;SLC25A3 PHC OK/SW-cl.48;Phosphate carrier protein;mitochondrial (Phosphate transport protein) (PTP) (Solute carrier family 25 member 3); mitochondrial (Phosphate transport protein) (PTP) (Solute carrier family 25 member 3)
Observed Band	
Calculated Molecular Weight	40kD
Cell Pathway	Mitochondrion inner membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

