



ZFPL1 mouse mAb

Catalog No	YP-mAb-10164
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZFPL1
Immunogen	Synthesized peptide derived from human ZFPL1 AA range: 240-290
Specificity	This antibody detects endogenous levels of human ZFPL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	zinc finger protein like 1;Zinc finger protein MCG4;Zinc finger protein-like 1; ZFPL1;Zinc finger protein-like 1 (Zinc finger protein MCG4)
Observed Band	35-40 kDa
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus, cis-Golgi network membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed strongly in the exocrine pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images