



ZNF41 mouse mAb

Catalog No	YP-mAb-10165
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF41
Immunogen	Synthesized peptide derived from human ZNF41 AA range: 667-717
Specificity	This antibody detects endogenous levels of human ZNF41
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZNF416;ZNF41;Zinc finger protein 41
Observed Band	67kDa
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in the heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images