



NICN1 mouse mAb

Catalog No	YP-mAb-10167
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NICN1
Immunogen	Synthesized peptide derived from human NICN1 AA range: 145-195
Specificity	This antibody detects endogenous levels of human NICN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NICN1;nicolin 1;PGs5;Nicolin-1 (NPCEDRG) (Tubulin polyglutamylase complex subunit 5) (PGs5)
Observed Band	27 kDa
Calculated Molecular Weight	24 kDa
Cell Pathway	Nucleus .
Tissue Specificity	High expression level is found in brain, testis, liver and kidney. Weak expression in spleen, leukocytes, small intestine and colon.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images