



MFS2B mouse mAb

Catalog No	YP-mAb-10169
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MFSD2B
Immunogen	Synthesized peptide derived from human MFS2B AA range: 145-195
Specificity	This antibody detects endogenous levels of human MFS2B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MFS2B;MFSD2B;Major facilitator superfamily domain-containing protein 2B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Localizes to the cell membrane and intracellular membranes. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images