



U1SBP mouse mAb

Catalog No	YP-mAb-10171
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SNRNP35 HM1 U1SNRNPBP
Immunogen	Synthesized peptide derived from human U1SBP AA range: 168-218
Specificity	This antibody detects endogenous levels of human U1SBP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	U1SBP;SNRNP35 HM1 U1SNRNPBP;U11/U12 small nuclear ribonucleoprotein 35 kDa protein (U11/U12 snRNP 35 kDa protein) (U11/U12-35K) (Protein HM-1) (U1 snRNP-binding protein homolog)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in heart, liver, skeletal muscle and pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images