



KIAS3 mouse mAb

Catalog No	YP-mAb-10178
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIRREL3-AS3 NCRNA00288 PRR10
Immunogen	Synthesized peptide derived from human KIAS3 AA range: 62-112
Specificity	This antibody detects endogenous levels of human KIAS3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAS3;KIRREL3-AS3 NCRNA00288 PRR10;Putative uncharacterized protein KIRREL3-AS3 (KIRREL3 antisense RNA 1) (KIRREL3 antisense gene protein 1) (Proline-rich protein 10)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images