



SUGT1 mouse mAb

Catalog No	YP-mAb-10183
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SUGT1
Immunogen	Synthesized peptide derived from human SUGT1 AA range: 96-146
Specificity	This antibody detects endogenous levels of human SUGT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Sgt1;Protein SGT1 homolog;Protein 40-6-3;Protein 40 6 3;SUGT1;Suppressor of G2 allele of SKP1 homolog (Protein 40-6-3) (Sgt1)
Observed Band	38 kDa and 41 kDa
Calculated Molecular Weight	41 kDa
Cell Pathway	Cytoplasm . Nucleus . Translocates to the nucleus upon heat shock, requiring S100A6.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images