



MOSC2 mouse mAb

Catalog No	YP-mAb-10184
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MARC2 MOSC2
Immunogen	Synthesized peptide derived from human MOSC2 AA range: 140-190
Specificity	This antibody detects endogenous levels of human MOSC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:1.7.-.-;mARC2;Mitochondrial amidoxime reducing component 2;Moco sulfurase C-terminal domain-containing protein 2;Molybdenum cofactor sulfurase C-terminal domain-containing protein 2;MOSC2;MARC2 MOSC2; MOSC domain-containing protein 2;mitochondrial (EC 1.-.-.-) (Mitochondrial amidoxime reducing component 2) (mARC2) (Moco sulfurase C-terminal domain-containing protein 2) (Molybdenum cofactor sulfura; mitochondrial (EC 1.-.-.-) (Mitochondrial amidoxime reducing component 2) (mARC2) (Moco sulfurase C-terminal domain-containing protein 2) (Molybdenum cofactor sulfura
Observed Band	35-38 kDa
Calculated Molecular Weight	335 aa, 38 kDa
Cell Pathway	Mitochondrion outer membrane; Peripheral membrane protein. Peroxisome .
Tissue Specificity	
Function	
Background	

**matters needing
attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images