



TM158 mouse mAb

Catalog No	YP-mAb-10189
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM158 HBBP RIS1
Immunogen	Synthesized peptide derived from human TM158 AA range: 97-147
Specificity	This antibody detects endogenous levels of human TM158
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	40 kDa BINP binding protein;40 kDa BINP-binding protein;BBP;BINP receptor; Brain injury derived neurotrophic peptide (BINP) binding protein;Brain specific binding protein;DKFZp586E1621;HBBP;p40BBP;Ras induced senescence protein 1;Ras-induced senescence protein 1;RIS1;TM158_HUMAN;TMEM 158; TMEM158;Transmembrane protein 158;TM158;TMEM158 HBBP RIS1; Transmembrane protein 158 (40 kDa BINP-binding protein) (p40BBP) (Ras-induced senescence protein 1)
Observed Band	30 kDa;
Calculated Molecular Weight	30kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images