



MARHB mouse mAb

Catalog No	YP-mAb-10193
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	44266
Immunogen	Synthesized peptide derived from human MARHB AA range: 198-248
Specificity	This antibody detects endogenous levels of human MARHB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MARHB;44266;E3 ubiquitin-protein ligase MARCH11 (EC 6.3.2.-) (Membrane-associated RING finger protein 11) (Membrane-associated RING-CH protein XI) (MARCH-XI)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images