



AP3B1 mouse mAb

Catalog No	YP-mAb-10194
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AP3B1 ADTB3A
Immunogen	Synthesized peptide derived from human AP3B1 AA range: 525-575
Specificity	This antibody detects endogenous levels of human AP3B1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADTB3;ADTB3A;AP 3 complex subunit beta 1;AP3B1;Beta 3A adaptin;HPS; HPS2;PE;AP3B1 ADTB3A;AP-3 complex subunit beta-1 (Adapter-related protein complex 3 subunit beta-1) (Adaptor protein complex AP-3 subunit beta-1) (Beta-3A-adaptin) (Clathrin assembly protein complex 3 beta-1 large chain)
Observed Band	140 kDa
Calculated Molecular Weight	1094 aa, 121 kDa
Cell Pathway	Cytoplasmic vesicle, clathrin-coated vesicle membrane ; Peripheral membrane protein ; Cytoplasmic side . Golgi apparatus . Component of the coat surrounding the cytoplasmic face of coated vesicles located at the Golgi complex. .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images