



CKP2L mouse mAb

Catalog No	YP-mAb-10197
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CKAP2L
Immunogen	Synthesized peptide derived from human CKP2L AA range: 349-399
Specificity	This antibody detects endogenous levels of human CKP2L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CKP2L;CKAP2L;Cytoskeleton-associated protein 2-like
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, spindle pole . Uniformly distributed along each microtubule bundle of spindles in addition to centrioles during mitosis, expression promptly diminishes at interphase. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images