



PRDM4 mouse mAb

Catalog No	YP-mAb-10201
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRDM4 PFM1
Immunogen	Synthesized peptide derived from human PRDM4 AA range: 85-135
Specificity	This antibody detects endogenous levels of human PRDM4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PFM1;PR domain containing 4;PR domain containing protein 4;PRDM4;PRDM4 PFM1;PR domain zinc finger protein 4 (PR domain-containing protein 4)
Observed Band	88 kDa
Calculated Molecular Weight	801 aa, 88 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in many tissues. Highly expressed in ovary, testis, pancreas, brain, heart and prostate.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images