



HEAT2 mouse mAb

Catalog No	YP-mAb-10202
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HEATR2
Immunogen	Synthesized peptide derived from human HEAT2 AA range: 242-292
Specificity	This antibody detects endogenous levels of human HEAT2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HEAT2;HEATR2;HEAT repeat-containing protein 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Dynein axonemal particle . Observed only in the cytoplasm of ciliated cells and absent from cilia. .
Tissue Specificity	Expressed in nasal epithelium and lung epithelium by ciliated cells (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images