



GID8 mouse mAb

Catalog No	YP-mAb-10203
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GID8 C20orf11
Immunogen	Synthesized peptide derived from human GID8 AA range: 115-165
Specificity	This antibody detects endogenous levels of human GID8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C20orf11;Glucose-induced degradation protein 8 homolog;Protein C20orf11;Twa1;Two hybrid-associated protein 1 with RanBPM;GID8;GID8 C20orf11;Glucose-induced degradation protein 8 homolog (Two hybrid-associated protein 1 with RanBPM) (Twa1)
Observed Band	30 kDa
Calculated Molecular Weight	228 aa, 27 kDa
Cell Pathway	Cytoplasm . Nucleus . Localizes in the cytoplasm in the absence of Wnt stimulation and in the nucleus in the presence of Wnt stimulation. .
Tissue Specificity	Up-regulated in colorectal cancer tissues (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images