



S29A3 mouse mAb

Catalog No	YP-mAb-10205
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC29A3 ENT3 UNQ717/PRO1380
Immunogen	Synthesized peptide derived from human S29A3
Specificity	This antibody detects endogenous levels of human S29A3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S29A3;SLC29A3 ENT3 UNQ717/PRO1380;Equilibrative nucleoside transporter 3 (hENT3) (Solute carrier family 29 member 3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane; Multi-pass membrane protein. Late endosome membrane. Lysosome membrane. Observed in a punctate intracellular pattern showing partial colocalization with late endosomes/lysosomes. Not detected at the cell surface.
Tissue Specificity	Widely expressed in both adult and fetal tissues. Highest levels in placenta, uterus, ovary, spleen, lymph node and bone marrow. Lowest levels in brain and heart.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images