



SDIM1 mouse mAb

Catalog No	YP-mAb-10208
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SDIM1
Immunogen	Synthesized peptide derived from human SDIM1 AA range: 83-133
Specificity	This antibody detects endogenous levels of human SDIM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SDIM1;Stress-responsive DNAJB4-interacting membrane protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in brain with higher detection in neurons than astrocytes. Decreased expression in Alzheimer brains. Detected at protein level in brain and cervix.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images