



AP4S1 mouse mAb

Catalog No	YP-mAb-10209
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AP4S1
Immunogen	Synthesized peptide derived from human AP4S1 AA range: 69-119
Specificity	This antibody detects endogenous levels of human AP4S1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AP4S1;AP-4 complex subunit sigma-1 (AP-4 adapter complex subunit sigma-1) (Adapter-related protein complex 4 subunit sigma-1) (Sigma-1 subunit of AP-4) (Sigma-4-adaptin) (Sigma4-adaptin)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images