



ANGE2 mouse mAb

Catalog No	YP-mAb-10211
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ANGEL2 KIAA0759L
Immunogen	Synthesized peptide derived from human ANGE2 AA range: 446-496
Specificity	This antibody detects endogenous levels of human ANGE2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANGE2_HUMAN;Angel homolog 2;angel2;FLJ12793;KIAA0759L;LOC90806 protein;OTTHUMP00000035020;OTTHUMP00000035021;Protein angel homolog 2;ANGE2;ANGEL2 KIAA0759L
Observed Band	62kD
Calculated Molecular Weight	62kD(
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images