



ECM29 mouse mAb

Catalog No	YP-mAb-10213
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ECM29 KIAA0368
Immunogen	Synthesized peptide derived from human ECM29 AA range: 1599-1649
Specificity	This antibody detects endogenous levels of human ECM29
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ECM29;extracellular mutant 29;KIAA0368;KIAA1962;RP11 386D8.2;ECM29 KIAA0368;Proteasome-associated protein ECM29 homolog (Ecm29)
Observed Band	205 kDa
Calculated Molecular Weight	1845 aa, 204 kDa
Cell Pathway	Endoplasmic reticulum . Endoplasmic reticulum-Golgi intermediate compartment . Endosome . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Nucleus . Endosome, multivesicular body . Cytoplasmic vesicle .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images