



SRP19 mouse mAb

Catalog No	YP-mAb-10214
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SRP19
Immunogen	Synthesized peptide derived from human SRP19 AA range: 22-72
Specificity	This antibody detects endogenous levels of human SRP19
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	signal recognition particle 19kDa;SRP 19;SRP19;Signal recognition particle 19 kDa protein (SRP19)
Observed Band	18-25 kDa
Calculated Molecular Weight	144 aa, 16 kDa
Cell Pathway	Cytoplasm . Nucleus, nucleolus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images