



PSMG1 mouse mAb

Catalog No	YP-mAb-10215
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PSMG1 C21LRP DSCR2 PAC1
Immunogen	Synthesized peptide derived from human PSMG1 AA range: 144-194
Specificity	This antibody detects endogenous levels of human PSMG1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PSMG1;PSMG1 C21LRP DSCR2 PAC1;Proteasome assembly chaperone 1 (PAC-1) (Chromosome 21 leucine-rich protein) (C21-LRP) (Down syndrome critical region protein 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm. Endoplasmic reticulum.
Tissue Specificity	In the adult, detected in brain, colon, leukocytes, breast and testis. Widely expressed in the fetus. Also expressed in a variety of proliferating cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images