



MRP4 mouse mAb

Catalog No	YP-mAb-10216
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABCC4 MRP4
Immunogen	Synthesized peptide derived from human MRP4 AA range: 884-934
Specificity	This antibody detects endogenous levels of human MRP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABCC4;EST170205;MOAT B;MOATB;MRP4;ABCC4 MRP4;Multidrug resistance-associated protein 4 (ATP-binding cassette sub-family C member 4) (MRP/cMOAT-related ABC transporter) (Multi-specific organic anion transporter B) (MOAT-B)
Observed Band	150-180 kDa
Calculated Molecular Weight	859 aa, 97 kDa
Cell Pathway	Basolateral cell membrane ; Multi-pass membrane protein . Apical cell membrane ; Multi-pass membrane protein . Its localization to the basolateral or apical membranes is tissue-dependent. .
Tissue Specificity	Widely expressed, with particularly high levels in prostate, but is barely detectable in liver. sinusoidal membrane of hepatocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images