



ZN451 mouse mAb

Catalog No	YP-mAb-10217
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZN451 COASTER KIAA0576 KIAA1702
Immunogen	Synthesized peptide derived from human ZN451 AA range: 222-272
Specificity	This antibody detects endogenous levels of human ZN451
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN451;ZNF451 COASTER KIAA0576 KIAA1702;Zinc finger protein 451 (Coactivator for steroid receptors)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Nucleus, PML body . Nucleus, nucleoplasm. Colocalizes with UBE2I/ UBC9, SUMO1 and SUMO2 in nuclear granules; this probably requires sumoylation. Desumoylation leads to diffuse nucleoplasmic location. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images