



SPTA3 mouse mAb

Catalog No	YP-mAb-10218
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPATA3 TSARG1
Immunogen	Synthesized peptide derived from human SPTA3 AA range: 37-87
Specificity	This antibody detects endogenous levels of human SPTA3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1700011N12Rik;1700029H01Rik;4930424D10Rik;Mtsarg1;SPATA3; Spermatogenesis-associated protein 3;SPTA3_HUMAN;Testis and spermatogenesis cell-related protein 1;Testis spermatocyte apoptosis-related protein 1;Tsarg1;SPTA3;SPATA3 TSARG1;Spermatogenesis-associated protein 3 (Testis and spermatogenesis cell-related protein 1) (Testis spermatocyte apoptosis-related protein 1)
Observed Band	21kD
Calculated Molecular Weight	21kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images