



NOMO2 mouse mAb

Catalog No	YP-mAb-10219
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NOMO2
Immunogen	Synthesized peptide derived from human NOMO2 AA range: 139-189
Specificity	This antibody detects endogenous levels of human NOMO2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	pM5 protein 2;Nodal modulator 2;BOS complex subunit NOMO2;NOMO2;Nodal modulator 2 (pM5 protein 2)
Observed Band	130 kDa
Calculated Molecular Weight	139 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Tissue Specificity	Highly expressed in pancreas and skeletal muscle and, at lower levels, in heart.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images