



CAR17 mouse mAb

Catalog No	YP-mAb-10220
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CARD17 INCA
Immunogen	Synthesized peptide derived from human CAR17 AA range: 22-72
Specificity	This antibody detects endogenous levels of human CAR17
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAR17;CARD17 INCA;Caspase recruitment domain-containing protein 17 (Caspase-1 inhibitor INCA) (Inhibitory caspase recruitment domain protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images