



ENK11 mouse mAb

Catalog No	YP-mAb-10222
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Synthesized peptide derived from human ENK11 AA range: 64-114
Specificity	This antibody detects endogenous levels of human ENK11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ENK11;HERV-K_22q11.21 provirus ancestral Env polyprotein (Envelope polyprotein) (HERV-K101 envelope protein) [Cleaved into: Surface protein (SU); Transmembrane protein (TM)]
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Transmembrane protein]: Cell membrane ; Single-pass type I membrane protein . ; [Surface protein]: Cell membrane ; Peripheral membrane protein . The surface protein is not anchored to the membrane, but localizes to the extracellular surface through its binding to TM. . ; [Endogenous retrovirus group K member 24 Env polyprotein]: Virion .
Tissue Specificity	Expressed in teratocarcinoma cell line.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images