



ISOC2 mouse mAb

Catalog No	YP-mAb-10225
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ISOC2
Immunogen	Synthesized peptide derived from human ISOC2 AA range: 47-97
Specificity	This antibody detects endogenous levels of human ISOC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ISOC2;Isochorismatase domain-containing protein 2;mitochondrial; mitochondrial
Observed Band	28 kDa
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Localizes to the nucleus in the presence of CDKN2A.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images