



CR032 mouse mAb

Catalog No	YP-mAb-10226
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C18orf32
Immunogen	Synthesized peptide derived from human CR032 AA range: 14-64
Specificity	This antibody detects endogenous levels of human CR032
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CR032;C18orf32;UPF0729 protein C18orf32 (Putative NF-kappa-B-activating protein 200)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum . Lipid droplet .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images