



DUS11 mouse mAb

Catalog No	YP-mAb-10228
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DUSP11 PIR1
Immunogen	Synthesized peptide derived from human DUS11 AA range: 218-268
Specificity	This antibody detects endogenous levels of human DUS11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Dual specificity phosphatase 11;Dual specificity protein phosphatase 11; DUS11_HUMAN;DUSP11;Phosphatase that interacts with RNA/RNP complex 1; PIR1;RNA/RNP complex 1 interacting;RNA/RNP complex interacting phosphatase;RNA/RNP complex-1-interacting phosphatase;RNA/RNP complex1 intereracting phosphatase;Serine/threonine specific protein phosphatase;DUS11;DUSP11 PIR1;RNA/RNP complex-1-interacting phosphatase (EC 3.1.3.-) (Dual specificity protein phosphatase 11) (Phosphatase that interacts with RNA/RNP complex 1)
Observed Band	39 kDa
Calculated Molecular Weight	44kD
Cell Pathway	Nucleus . Nucleus speckle .
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images