



EIF3M mouse mAb

Catalog No	YP-mAb-10229
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EIF3M HFLB5 PCID1 GA17 PNAS-125
Immunogen	Synthesized peptide derived from human EIF3M AA range: 107-157
Specificity	This antibody detects endogenous levels of human EIF3M
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PCID1;hFL-B5;HFLB5;hfl B5;GA17;EIF3M;EIF3M HFLB5 PCID1 GA17 PNAS-125;Eukaryotic translation initiation factor 3 subunit M (eIF3m) (Fetal lung protein B5) (hFL-B5) (PCI domain-containing protein 1)
Observed Band	35-43 kDa
Calculated Molecular Weight	374 aa, 43 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Broadly expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images