



PAQRB mouse mAb

Catalog No	YP-mAb-10230
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MMD PAQR11
Immunogen	Synthesized peptide derived from human PAQRB AA range: 38-88
Specificity	This antibody detects endogenous levels of human PAQRB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAQRB;MMD PAQR11;Monocyte to macrophage differentiation protein (Progestin and adipoQ receptor family member 11) (Progestin and adipoQ receptor family member XI)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Late endosome membrane ; Multi-pass membrane protein . Lysosome membrane ; Multi-pass membrane protein .
Tissue Specificity	Exhibits relatively ubiquitous expression with preferential expression in mature (in vitro differentiated) macrophages.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images