



PANK3 mouse mAb

Catalog No	YP-mAb-10231
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PANK3
Immunogen	Synthesized peptide derived from human PANK3 AA range: 82-132
Specificity	This antibody detects endogenous levels of human PANK3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hPank3;PANK1/2/3;PANK3;pantothenate kinase 3;Pantothenic acid kinase 3; Pantothenate kinase 3 (hPank3) (EC 2.7.1.33) (Pantothenic acid kinase 3)
Observed Band	41 kDa
Calculated Molecular Weight	370 aa, 41 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Highly expressed in the liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images