



TTL11 mouse mAb

Catalog No	YP-mAb-10232
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TTLL11 C9orf20
Immunogen	Synthesized peptide derived from human TTL11 AA range: 158-208
Specificity	This antibody detects endogenous levels of human TTL11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TTL11;TTLL11 C9orf20;Tubulin polyglutamylase TTLL11 (EC 6.-.-.) (Tubulin--tyrosine ligase-like protein 11)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, cilium basal body . Cytoplasm, cytoskeleton .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images