



IFM5 mouse mAb

Catalog No	YP-mAb-10233
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IFITM5
Immunogen	Synthesized peptide derived from human IFM5 AA range: 1-51
Specificity	This antibody detects endogenous levels of human IFM5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IFM5;IFITM5;Interferon-induced transmembrane protein 5 (Bone-restricted interferon-induced transmembrane protein-like protein) (BRIL) (Dispanin subfamily A member 1) (DSPA1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Detected in bone (PubMed:24058703). Detected in osteoblasts and fibroblasts (at protein level) (PubMed:24519609). Detected in bone (PubMed:24058703). Detected in osteoblasts and fibroblasts (PubMed:24519609).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images