



S15A5 mouse mAb

Catalog No	YP-mAb-10234
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC15A5
Immunogen	Synthesized peptide derived from human S15A5 AA range: 430-480
Specificity	This antibody detects endogenous levels of human S15A5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S15A5;SLC15A5;Solute carrier family 15 member 5
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images