



EGFEM mouse mAb

Catalog No	YP-mAb-10235
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EGFEM1P C3orf50 NCRNA00259
Immunogen	Synthesized peptide derived from human EGFEM AA range: 33-83
Specificity	This antibody detects endogenous levels of human EGFEM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EGFEM;EGFEM1P C3orf50 NCRNA00259;Putative EGF-like and EMI domain-containing protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images