



DNM3L mouse mAb

Catalog No	YP-mAb-10239
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNMT3L
Immunogen	Synthesized peptide derived from human DNM3L AA range: 52-102
Specificity	This antibody detects endogenous levels of human DNM3L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DNM3L;DNMT3L;DNA (cytosine-5)-methyltransferase 3-like
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed at low levels in several tissues including testis, ovary, and thymus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

