



TNPO2 mouse mAb

Catalog No	YP-mAb-10240
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TNPO2
Immunogen	Synthesized peptide derived from human TNPO2 AA range: 306-356
Specificity	This antibody detects endogenous levels of human TNPO2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TNPO 2;KPNB2B;Karyopherin beta-2b;IPO3;Importin 3;TNPO2;Transportin-2 (Karyopherin beta-2b)
Observed Band	90-100 kDa
Calculated Molecular Weight	887 aa, 100 kDa
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images