



F108C mouse mAb

Catalog No	YP-mAb-10241
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM108C1
Immunogen	Synthesized peptide derived from human F108C AA range: 233-283
Specificity	This antibody detects endogenous levels of human F108C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F108C;FAM108C1;Abhydrolase domain-containing protein FAM108C1 (EC 3.-.-)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Recycling endosome membrane ; Lipid-anchor ; Cytoplasmic side . Cell projection, dendritic spine . Cell junction, synapse, postsynaptic density membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images