



ZKSC3 mouse mAb

Catalog No	YP-mAb-10242
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZKSCAN3 ZFP47 ZNF306 ZNF309 ZSCAN13
Immunogen	Synthesized peptide derived from human ZKSC3 AA range: 457-507
Specificity	This antibody detects endogenous levels of human ZKSC3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZKSC3;ZKSCAN3 ZFP47 ZNF306 ZNF309 ZSCAN13;Zinc finger protein with KRAB and SCAN domains 3 (Zinc finger and SCAN domain-containing protein 13) (Zinc finger protein 306) (Zinc finger protein 309) (Zinc finger protein 47 homolog) (Zf47) (Zfp-47
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus. Cytoplasm. Mainly localizes in the nucleus. Under starvation conditions translocates to the cytoplasm, allowing expression of target genes involved in autophagy and lysosome biogenesis/function.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images